

Product Guide

Mark II™

Trusted by many for reliable
performance and excellent quality



Features & benefits:

- Excellent atomisation capabilities with great performance and quality
- Has a Drug Master File (DMF)
- Optional configuration for upside-down spraying
- Standard white colour; other colours available upon request
- Available with shipping clip for various sizes and outputs
- Available with protective hood for nasal adapters in custom colours

Adapters

Nasal,
Oral,
Topical,
Ear



Round



Square



Oval²



Mini



Elegant



Shroud




Hexagon



Throat adapter short



Throat adapter medium




Throat adapter long



Throat adapter long with wings



Throat adapter articulated¹



ATA 19-H



Dispensing head



Spray head



Spray head with capillary tube²



Ear adapter

Closures



Screw
(GL18 short, 18-415,
GL20, 18-410,
22-415, GL18,
PP18S)



Ferrule
(20 mm-gold
or silver)



Snap
(S20, GS20)

Overcaps



Nasal hoods



Throat adapter
medium
protection

Options



Safety clip²

Output

0.100 ml, 0.120 ml,
0.150 ml, 0.160 ml,
0.180 ml, 0.200 ml

Cross-sectional drawing for Mark II is available at SilganDS.com

¹ compatible with ferrule closures

² available upon special request